

NEW
from Matrix Separations

Introducing the Matrix HD-Sb Hydrogel Membrane

*Cation Exchanger augmented with HIC Groups
for Purification of Monoclonal Antibodies*

- Extremely fast Capture or Polishing
- Very high dynamic binding capacity, excellent monomer recovery and high aggregate clearance
- High salt tolerance for robust processes



Explore the future of disposable downstream processing!

The new Recon HD-Sb membrane device is now available for lab scale exploration and DOE work. Low pressure Recon devices, ready-to-use with Luer lock fittings, are compatible with all benchtop chromatography systems. Larger devices for process development and clinical manufacturing are currently under development.

High Binding Capacity and high speed polishing with proven Matrix HD hydrogel technology

The high ligand density and excellent mass transfer of the three-dimensional macroporous hydrogel give a binding capacity >90 mg mAb/ml of membrane (10% breakthrough) at flow rates of 10 MV/min or greater (residence time of 6 seconds or less). In flowthrough mode, a single 0.8 ml Recon device can polish 250 mg of monoclonal antibody in one step.

High Aggregate Clearance and HCP Removal by mixed-mode chemistries

The HD-Sb membrane is a mixture of strong cation exchange (sulfonic acid) and hydrophobic (t-butyl) groups. This unique chemistry achieves highly purified antibodies in both bind & elute and flowthrough modes.

LEARN MORE: For more information, please visit: www.natrixseparations.com/Sb




www.JSRLifeSciences.com

Distributed globally by JSR Life Sciences

JSR Life Sciences 
MATERIALS INNOVATION

Natrix HD-Sb Hydrogel Membrane Columns: Performance Data



Figures A & B: High Aggregate Removal in Bind & Elute and Flowthrough Mode

SEC analyses of Feed and Eluate from Bind & Elute [Figure A] and Feed and Flowthrough [Figure B] experiments show highly efficient removal of aggregates. Load (50 mg/mL [A], 300 mg/mL [B]), flow rate 10 membrane volumes/minute, sample concentration (1 g/L [A], 3 g/L [B]) in equilibration buffer (45 mM Na-acetate, 130 mM NaCl, pH4.5 [A] pH5 [B], 16.3 mS/cm [A] 10 mS/cm [B]). Elution with wash buffer (20 mM phosphate, pH6.3) at 6 mS/cm.

Table 1: High Monomer yield in Bind & Elute and Flowthrough modes

The table below shows the exceptional clearance of aggregates and high monomer yield in both Bind & Elute and Flowthrough mode when used for polishing antibodies purified using Protein A chromatography.

| Mode | Load | Aggregates | | Monomer Yield |
|--------------|-----------|------------|--------------|---------------|
| | | Feed | Elution | |
| Bind & Elute | 50 mg/mL | 1.8% | Undetectable | >90% |
| Bind & Elute | 50 mg/mL | 3.3% | Undetectable | >90% |
| Bind & Elute | 50 mg/mL | 8.1% | Undetectable | >90% |
| | | Feed | Flowthrough | |
| Flowthrough | 300 mg/mL | 12.8% | 0.6% | >90% |

Same antibody used in all experiments. Different Aggregates levels created by changing the time of the low pH holding step.

Ordering and Customer Service Support

Please contact our global distributor, JSR Life Sciences:

| REGION | | |
|-----------------|--------------------|------------------|
| North America | Asia | Europe |
| cs@jsrmicro.com | info@jls.jsr.co.jp | bio@jsrmicro.be |
| +1.408.543.8869 | +81.29.856.1106 | +32 (16) 832 832 |

Technical Inquiries:

Contact Natrix Separations:

+1 905.319.2682

info@natrixseparations.com

JULY 2015

Table 2: High HCP and aggregate reduction in capture and intermediate applications.

The HD-Sb is capable of high impurities clearance for purification of monoclonal antibodies. Two different monoclonal antibodies were used for capture studies and for bind & elute studies after a ProA affinity purification step. The HD-Sb shows good clearance for HCP and aggregates in both cases.

| Mode | Capture | Intermediate (B/E) |
|---------------------|--------------|--------------------|
| Load (mg/mL) | 55 | 65 |
| Flow rate (MV/min) | 10 | |
| HCP Clearance (LRV) | 1.4 | 1.7 |
| Load (ppm) | 350,800 | 140 |
| Eluate (ppm) | 12,730 | 2.8 |
| Aggregate Clearance | | |
| Load | Not measured | 1.80% |
| Eluate | 2.60% | Not detectable |
| mAb Yield | > 90% | > 90% |



Natrix Separations, Inc.

5295 John Lucas Drive | Unit 6

Burlington, Ontario L7L 6A8 CANADA

+1 905.319.2682 (Direct)